

Recombinant Human VEGF121

Cat. #: SVF1-010

Product Specifications

- Expression of Human Proteins in Human Cells
- Extremely Low Endotoxin Level
- High Purity
- Animal Free and Xeno Free
- Tag Free

Protein Description

Vascular endothelial growth factor (VEGF) is a potent growth and angiogenic cytokine. It is a member of the PDGF family that is characterized by the presence of eight conserved cysteine residues and a cystine knot structure. VEGF stimulates angiogenesis, vasculogenesis and endothelial cell growth, induces endothelial cell proliferation, promotes cell migration, inhibits apoptosis and induces permeabilization of blood vessels. Human express alternately spliced isoforms, and VEGF121 is the only form that lacks a basic heparin binding region and is freely diffusible. VEGF121 induces the proliferation of lymphatic endothelial cells and the lymph angiogenesis may be promoted by up regulation of VEGF121.

References

Leung DW, et al. (1989) Science 246, 1306-1309.

Robinson CJ, et al. (2001) J. Cell. Sci. 114, 853.

Tischer E, et al. (1991) J. Biol. Chem. 266, 11947-11954.

Byrne AM, et al. (2005) J. Cell. Mol. Med. 9, 777.

Source: Derived from human cells

Size: 10 µg

Shipping: Ambient temperature

Structure: Glycosylated homodimer

Purity: >95% by SDS-PAGE

Endotoxin Level: <0.5 EU/µg

Molecular Weight: 37 kDa homodimer, 50 kDa homotrimer

Formulation: Lyophilized from a 0.2 µm filtered solution in PBS without carrier protein

Activity Assay

Activity was measured by its ability to stimulate the proliferation of HUVEC (human umbilical vein endothelial cells).

Reconstitution

Briefly centrifuge the vial before opening. It is recommended to reconstitute the protein in sterile PBS containing 0.1% endotoxin-free recombinant human serum albumin.

Stability and Storage

Use a manual defrost freezer and avoid repeated freeze-thaw cycles. In general: 12 months from date of receipt, -20 to -80° C as supplied. 1 month, 2 to 8° C under sterile conditions after reconstitution. 3 months, -20 to -80° C under sterile conditions after reconstitution.



2600 Hilltop Drive, Building B, Suite C328, Richmond, CA 94806

Tel: (510) 708-0096 | Fax: (866) 605-8766 | www.alstembio.com | info@alstembio.com